

**ULRICH**

*100 Jahre*  
den Menschen  
im Fokus

STERILE CONTAINER SYSTEM.

**ULRICH SWISS.**



## STERILE CONTAINER SYSTEM.

### ULRICH SWISS.

#### Use and benefit

From keeping surgical instruments in an optimal sterile condition during transport and storage, to their use in surgical applications. In addition to safety aspects, the **ULRICH sterile containers** convince from an economic perspective - thanks to their very positive price/performance ratio.

By using these sterile containers, the entire sterile goods cycle can be standardized, controlled and documented from supply to disposal. This applies to the sterile storage and transport of the instruments, as well as to their ease of provision in the operating theatre and their safe disposal after use. When disposing of instruments, the **ULRICH sterile containers** reliably protect the hospital staff against contamination, thereby supporting the hospital's quality management in terms of hygiene and infection prevention and control.

#### Handling

- The ULRICH Bio Barrier containers are easy and quick to pack
- The filterless system means a filter change doesn't get forgotten or avoids those lengthy checks and documentation of Teflon filters with max. 1,000 cycles
- The ULRICH Bio Barrier containers are very robust and durable

## Organization

- The option to use colour-coding (lid and labels) simplifies organisation
- The lids and container bottoms have a UDI code. This means the number of reprocessing cycles can be tracked with an instrument management system. These key figures support your life cycle management (life cycle costing / maintenance checks, etc.) and therefore offer you improved safety
- Silicone-coated handles (heat and noise protection) mean the instruments can be used shortly after sterilization

## Cost-effectiveness

- The ULRICH sterile containers have very positive price/performance ratio
- The containers don't have to be purchased, they can also be loaned for a defined period

## Ecology


- The ULRICH Bio Barrier sterile containers offer an ecological advantage in that the reusable sterile good packaging systems do not need any paper or Teflon filters over their entire life cycle

## Maintenance

- The ULRICH sterile containers come with a two-year warranty
- The sterile containers can be sent to ULRICH AG for repairs, maintenance or functional checks

Learn more [www.ulrich-swiss.ch](http://www.ulrich-swiss.ch)

**1/1 Model Bio-Barrier,**  
Lid non perforated, Bottom non perforated,  
Bottom aluminum, Lid anodized aluminum



**Lids**

*Length x width*  
580 x 280 mm

U95-1000-01	gray
U95-1000-02	yellow
U95-1000-03	green
U95-1000-04	blue
U95-1000-05	red
U95-1000-06	black

**Bottom**

*Outside dimensions Bottom*  
*Length x width x height*

U95-1100-00	600 x 272 x 100 mm
U95-1116-00	600 x 272 x 116 mm
U95-1138-00	600 x 272 x 138 mm
U95-1190-00	600 x 272 x 190 mm
U95-1235-00	600 x 272 x 235 mm

**1/2 Model Bio-Barrier,**

Lid non perforated, Bottom non perforated,  
Bottom aluminum, Lid anodized aluminum



**Lids**

*Length x width*  
308 x 280 mm

U95-2000-01	gray
U95-2000-02	yellow
U95-2000-03	green
U95-2000-04	blue
U95-2000-05	red
U95-2000-06	black

**Bottom**

*Outside dimensions Bottom*  
*Length x width x height*

U95-2100-00	305 x 272 x 100 mm
U95-2116-00	305 x 272 x 116 mm
U95-2138-00	305 x 272 x 138 mm
U95-2190-00	305 x 272 x 190 mm
U95-2235-00	305 x 272 x 235 mm

**3/4 Model Bio-Barrier,**

Lid non perforated, Bottom non perforated,  
Bottom aluminum, Lid anodized aluminum



**Lids**

*Length x width*  
465 x 280 mm

U95-3000-01	gray
U95-3000-02	yellow
U95-3000-03	green
U95-3000-04	blue
U95-3000-05	red
U95-3000-06	black

**Bottom**

*Outside dimensions Bottom*  
*Length x width x height*

U95-3100-00	485 x 272 x 100 mm
U95-3116-00	485 x 272 x 116 mm
U95-3138-00	485 x 272 x 138 mm

**Mini-Container**

Lid perforated, Bottom non perforated,  
Bottom aluminum, Lid anodized aluminum



**Lids**

*Length x width*  
320 x 140 mm

U95-4000-01	gray
U95-4000-02	yellow
U95-4000-03	green
U95-4000-04	blue
U95-4000-05	red
U95-4000-06	black

**Bottom**

*Outside dimensions Bottom*  
*Length x width x height*

U95-4032-00	315 x 132 x 32 mm
U95-4063-00	315 x 132 x 63 mm
U95-4092-00	315 x 132 x 92 mm

**Dental-Container**

Lid perforated, Bottom non perforated,  
Bottom aluminum, Lid anodized aluminum



**Lids**

*Length x width*  
325 x 192 mm

U95-4001-01	gray
U95-4001-02	yellow
U95-4001-03	green
U95-4001-04	blue
U95-4001-05	red
U95-4001-06	black

**Bottom**

*Outside dimensions Bottom*  
*Length x width x height*

U95-4043-00	320 x 186 x 43 mm
U95-4059-00	320 x 186 x 59 mm
U95-4122-00	320 x 186 x 122 mm



Additional models on request.

### Product range

The ULRICH sterile containers are available in various designs:

- Bio Barrier model
- Window model
- Conventional filter systems (paper or teflon)

All these container models come in standard sizes 1/1, 1/2, 3/4 and with different lid colours.

### ULRICH AG. ST. GALLEN.

MÖVENSTRASSE 12  
CH-9015 ST. GALLEN

T +41 71 314 62 62  
F +41 71 314 62 99

[www.ulrich-swiss.ch](http://www.ulrich-swiss.ch)  
[info@ulrich-swiss.ch](mailto:info@ulrich-swiss.ch)